# Package 'rxylib'

January 4, 2023
Type Package
Title Import XY-Data into R
<b>Description</b> Provides access to the 'xylib' C library for to import xy data from powder diffraction, spectroscopy and other experimental methods.
Version 0.2.9
<b>Date</b> 2023-01-04
Author Sebastian Kreutzer [aut, trl, cre] ( <a href="https://orcid.org/0000-0002-0734-2199">https://orcid.org/0000-0002-0734-2199</a> ), Johannes Friedrich [aut] ( <a href="https://orcid.org/0000-0002-0805-9547">https://orcid.org/0000-0002-0805-9547</a> ), RLum Team [ctb], Marcin Wojdyr [cph] (C++ library 'xylib'), Peng Zhang [cph] (C++ library 'xylib')
Maintainer Sebastian Kreutzer < sebastian.kreutzer@uni-heidelberg.de>
<pre>URL https://github.com/R-Lum/rxylib</pre>
<pre>BugReports https://github.com/R-Lum/rxylib/issues</pre>
License GPL-3   LGPL-2.1
<b>Depends</b> R ( $>= 4.1$ ), utils
Imports methods, Rcpp (>= 1.0.9)
Suggests testthat (>= 3.1.4)
<b>LinkingTo</b> Rcpp (>= 1.0.9), BH (>= 1.78.0)
Encoding UTF-8
Language en-GB
Collate 'methods_rxylib.R' 'rxylib.R' 'RcppExports.R' 'read_xyData.R' 'convert_xy2TKA.R'
RoxygenNote 7.2.3
NeedsCompilation yes
R topics documented:
rxylib-package convert_xy2TKA methods_rxylib read_xyData
Index

2 rxylib-package

# Description

Provides access to the 'xylib' C++ library for to import xy data from powder diffraction, spectroscopy and other experimental methods, like gamma-ray spectrometry.

License: GPL-3 | LGPL-2.1 (for the C++ library 'xylib')

# **Details**

# **Funding**

Between 2017-2019, the work of Sebastian Kreutzer as maintainer of the package was supported by LabEx LaScArBx (ANR - n. ANR-10-LABX-52).

From 01/2020-04/2022, Sebastian Kreutzer received funding from the European Union's Horizon 2020 research and innovation programme under the Marie Skłodowska-Curie grant agreement No 844457 (project: CREDit).

Supported data formats: library version: 1.6.0

IDNAMEDESCRIPTIONFILE EXTENSIONVALID_OPTIONSDATATYPI[1,]cpiSietronics Sieray CPIcpiascii[2,]uxdBruker Diffrac-AT UXDuxdascii[3,]rigaku_datRigaku DATdatascii[4,]bruker_rawSiemens/Bruker RAWrawbinary[5,]bruker_spcBruker ESP300-E SPCspcbinary[6,]vamasVAMAS ISO-14976vmsascii[7,]philips_udfPhilips UDFudfascii[8,]spePI WinSpec SPEspebinary[9,]pdcifPowder Diffraction CIFcifascii[10,]philips_rdPhilips PC-APD RD/SDrd sdbinary[11,]xrdmlPANalytical XRDMLxrdmlascii[12,]canberra_mcaCanberra MCAmca datbinary	E <b>BI</b> sin mı mı
[2,] uxd Bruker Diffrac-AT UXD uxd ascii [3,] rigaku_dat Rigaku DAT dat ascii [4,] bruker_raw Siemens/Bruker RAW raw binary [5,] bruker_spc Bruker ESP300-E SPC spc binary [6,] vamas VAMAS ISO-14976 vms ascii [7,] philips_udf Philips UDF udf ascii [8,] spe PI WinSpec SPE spe binary [9,] pdcif Powder Diffraction CIF cif ascii [10,] philips_rd Philips PC-APD RD/SD rd sd binary [11,] xrdml PANalytical XRDML xrdml ascii	mι
[3,] rigaku_dat Rigaku DAT dat ascii [4,] bruker_raw Siemens/Bruker RAW raw binary [5,] bruker_spc Bruker ESP300-E SPC spc binary [6,] vamas VAMAS ISO-14976 vms ascii [7,] philips_udf Philips UDF udf ascii [8,] spe PI WinSpec SPE spe binary [9,] pdcif Powder Diffraction CIF cif ascii [10,] philips_rd Philips PC-APD RD/SD rd sd binary [11,] xrdml PANalytical XRDML xrdml ascii	
[4,] bruker_raw Siemens/Bruker RAW raw binary [5,] bruker_spc Bruker ESP300-E SPC spc binary [6,] vamas VAMAS ISO-14976 vms ascii [7,] philips_udf Philips UDF udf ascii [8,] spe PI WinSpec SPE spe binary [9,] pdcif Powder Diffraction CIF cif ascii [10,] philips_rd Philips PC-APD RD/SD rd sd binary [11,] xrdml PANalytical XRDML xrdml ascii	mı
[5,]bruker_spcBruker ESP300-E SPCspcbinary[6,]vamasVAMAS ISO-14976vmsascii[7,]philips_udfPhilips UDFudfascii[8,]spePI WinSpec SPEspebinary[9,]pdcifPowder Diffraction CIFcifascii[10,]philips_rdPhilips PC-APD RD/SDrd sdbinary[11,]xrdmlPANalytical XRDMLxrdmlascii	
[6,]vamasVAMAS ISO-14976vmsascii[7,]philips_udfPhilips UDFudfascii[8,]spePI WinSpec SPEspebinary[9,]pdcifPowder Diffraction CIFcifascii[10,]philips_rdPhilips PC-APD RD/SDrd sdbinary[11,]xrdmlPANalytical XRDMLxrdmlascii	mı
[7,]philips_udfPhilips UDFudfascii[8,]spePI WinSpec SPEspebinary[9,]pdcifPowder Diffraction CIFcifascii[10,]philips_rdPhilips PC-APD RD/SDrd sdbinary[11,]xrdmlPANalytical XRDMLxrdmlascii	sin
[8,]spePI WinSpec SPEspebinary[9,]pdcifPowder Diffraction CIFcifascii[10,]philips_rdPhilips PC-APD RD/SDrd sdbinary[11,]xrdmlPANalytical XRDMLxrdmlascii	mı
[9,]pdcifPowder Diffraction CIFcifascii[10,]philips_rdPhilips PC-APD RD/SDrd sdbinary[11,]xrdmlPANalytical XRDMLxrdmlascii	sin
[10,] philips_rd Philips PC-APD RD/SD rd sd binary [11,] xrdml PANalytical XRDML xrdml ascii	mι
[11,] xrdml PANalytical XRDML xrdml ascii	mι
· · · · · · · · · · · · · · · · · · ·	sin
[12,] canberra mca Canberra MCA mca dat binary	mι
	sin
[13,] canberra_cnf Canberra CNF cnf binary	sin
[14,] xfit_xdd XFIT XDD xdd ascii	sin
[15,] riet7 RIET7/LHPM/PSI_DMC dat ascii	sin
[16,] dbws DBWS data dbw rit neu ascii	sin
[17,] chiplot ChiPLOT data chi ascii	sin
[18,] spectra Spectra / VGX 900 1 2 3 4 5 6 7 8 9 ascii	mι
[19,] specsxy SPECS SpecsLab2 xy xy ascii	mu
[20,] csv CSV or TSV csv tsv tab decimal-comma ascii	sin
[21,] xsyg Freiberg Instruments XSYG xsyg ascii	mι

convert\_xy2TKA 3

#### Author(s)

Sebastian Kreutzer, Institute of Geography, Ruprecht-Karl-University of Heidelberg (Germany), Johannes Friedrich (University of Bayreuth, Germany), RLum Team (family support), Marcin Wojdyr (maintainer and author of the C++ library xylib), Peng Zhang (author of the C++ library xylib)

convert\_xy2TKA

Convert xy-data to TKA

#### **Description**

Convert data to the Toolkit file format (TKA) as exported by, e.g., by the software Canberra Genie 2000.

#### Usage

```
convert_xy2TKA(object, file = NULL, overwrite = FALSE)
```

# **Arguments**

object rxylib (required): xy data as imported by the function read\_xyData. Optional a

file supported by the rxylib-package can be provided as input. Arguments can

be provided as list.

file character (optional): optional file path or file name for the output to be written.

If only a path is provided the output file name is derived from the input file name.

Argument can be provided as list.

overwrite logical (with default): force overwriting of existing files if TRUE.

# **Details**

#### **Supported formats**

- · Canberra CNF
- further formats on request ...

#### Value

Returns a list of matrix objects or an output TKA-file.

#### **Function version**

0.1.1

# How to cite

Kreutzer, S., 2023. convert\_xy2TKA(): Convert xy-data to TKA. Function version 0.1.1. In: Kreutzer, S., Friedrich, J., 2023. rxylib: Import XY-Data into R . R package version 0.2.9. https://github.com/R-Lum/rxylib

#### Author(s)

Sebastian Kreutzer, Institute of Geography, Universität Heidelberg, Germany

4 methods\_rxylib

# **Examples**

```
##convert CNF data (no export to file system)
convert_xy2TKA(
  object = system.file("extdata/ExampleSpectrum.CNF", package = "rxylib"))

## Not run:
##export as file

##create temporary filepath
##(for usage replace by own path)
temp_file <- tempfile(pattern = "output", fileext = ".TKA")

##convert and write to file system
convert_xy2TKA(
  object = system.file("extdata/ExampleSpectrum.CNF", package = "rxylib"),
  file = temp_file)

## End(Not run)</pre>
```

methods\_rxylib

methods\_rxylib

# **Description**

S3-methods support by the package rxylib. Listed functions can be passed directly into S3 generics (e.g., plot, print) without reshaping the data.

# Usage

```
## S3 method for class 'rxylib'
print(x, ...)
## S3 method for class 'rxylib'
plot(x, block = NULL, ...)
```

# **Arguments**

```
    x (required): input object
    ... further arguments that can be passed to the method
    block numeric (with default): select block for plotting, e.g. c(1:2).
```

read\_xyData 5

read_xyData Import xy-Data for Supported Formats into R	read_xyData	Import xy-Data for Supported Formats into R	
---	-------------	---	--

# **Description**

The function provides access to the underlying xylib to import data for supported file formats into R. In most cases, only the file path is needed with further arguments to import the data. The function automatically recognises allowed formats. See rxylib-package for supported formats.

#### Usage

```
read_xyData(file, options = "", verbose = TRUE, metaData = TRUE)
```

#### **Arguments**

file	$\mbox{\it character}$ ( $\mbox{\it required}$ ): path and file to be imported. The argument accepts an URL.
options	character (with default): set format options (see rxylib-package)
verbose	logical (with default): enables/disables verbose mode
metaData	logical (with default): enables/disables the export of metadata

#### Value

The functions returns a list of matrices.

#### **Function version**

0.3.0

# How to cite

Kreutzer, S., Friedrich, J., 2023. read\_xyData(): Import xy-Data for Supported Formats into R. Function version 0.3.0. In: Kreutzer, S., Friedrich, J., 2023. rxylib: Import XY-Data into R. R package version 0.2.9. https://github.com/R-Lum/rxylib

# Author(s)

Sebastian Kreutzer, Institute of Geography, Ruprecht-Karl-University of Heidelberg (Germany), Johannes Friedrich, University of Bayreuth (Germany)

# **Examples**

```
##load example dataset
file <- system.file("extdata/ExampleSpectrum.CNF", package = "rxylib")
results <- read_xyData(file)
results

##plot xy-spectrum
plot(results,
    type = "1",
    xlab = "Energy [keV]",
    ylab = "Counts",</pre>
```

6 read\_xyData

```
main = "Thorite - 1800 s")

mtext(side = 3, "Canberra Inspector 1000, 3 x 3 NaI probe")

##plot contour for TL-spectrum

##imported from an XSYG-file
spectrum <- read_xyData(system.file("extdata/TLSpectrum.xsyg", package = "rxylib"))
contour(
    x = spectrum$dataset[[1]]$data_block[,1],
    y = 1:ncol(spectrum$dataset[[1]]$data_block[,-1]),
    z = spectrum$dataset[[1]]$data_block[,-1],
    xlab = "Wavelength [nm]",
    ylab = "#Channel",
    main = "TL Spectrum")</pre>
```

# **Index**

```
* IO
    convert_xy2TKA, 3
    read_xyData, 5
* package
    rxylib-package, 2
character, 3, 5
convert\_xy2TKA, 3
list, 3, 5
logical, 3, 5
matrix, 3
{\tt methods\_rxylib}, {\tt 4}
numeric, 4
plot, 4
\verb|plot.rxylib| (\verb|methods_rxylib|), 4
print,4
print.rxylib (methods_rxylib), 4
read_xyData, 3, 5
rxylib, 3
rxylib (rxylib-package), 2
rxylib-package, 2, 5
```